

HERE AT DR PHYSIO USA, WE MAKE THE MOST IMAGINATIVE, UNWINDING, AND RESTORATIVE, WELLBEING AND INDIVIDUAL CARE PRODUCTS.

In the hubbub of today's lifestyle, it becomes impossible to nurture yourself at times. In the chaos of maintaining every aspect of life efficiently, we forget that our bodies are not machines. They need absolute care and tenderness. Pain in neck, shoulder, back or foot should never be taken for granted. They can lead to serious ailments in the long run. To help you ease the situations of painfulness, Dr Physio USA brings to you excellent ways of relaxation and pain relief through our premium gadgets. We offer a non-intrusive answer to help renew and reestablish your muscles and relieve throbbing pain with our various products.

The variety of items and massagers we offer are intended to target particular areas of the body viz neck, shoulder, back, foot and so on. The adaptable capabilities of our gadgets are meant to suit your individual needs. The products we offer are meant for utmost body relaxation after a long day at work. Since we are continually in pursuit of excellence, our research and development team at NY, USA works diligently to ensure that each of our innovations are available at a value you would appreciate.

SIT BACK, RELAX AND ENJOY THE AMAZING DR PHYSIO EXPERIENCE!!



CLEOPATRA FACE CLEANING BRUSH



**USER
INSTRUCTIONS**

Thank you for purchasing Dr Physio Cleopatra Facial Cleaning Brush! Please read the instructions manual carefully before using it.

CONTENTS

1.Introduction	2
2.Key Attachments and Their Use	3
3.Setup Guide	5
4.FAQs	5
5.Benefits Of Using The Facial Brush	7
6.Using The Unit Safely	7
7.Warnings	8
8.Package Contents	9
9.Specifications	10
10.Customer Support	11

INTRODUCTION

With Dr Physio Cleopatra Facial Cleaning Brush, you can enjoy a comfortable yet powerful tapping massage therapy at home. It comes with 4 exchangeable attachments: big brush, small brush, pumice and face sponge which actively removes dirt, oil and dead cells from all over your face with deep cleansing of every pore. It works gently on all types of skin and reaches all those areas of your face where your hands cannot reach. It is very simple to use unit that helps in getting you healthy, squeaky-clean and smooth skin every time along with increasing the blood circulation.

Please read our instructions manual for the right use of Dr Physio Cleopatra Facial Cleaning Brush before starting. If you have any questions regarding any of the instructions or warnings, please contact our customer care service for help and support.

KEY ATTACHMENTS AND THEIR USE



Big Brush



Small Brush



Pumice Stone



Face Sponge



BIG BRUSH



Big brush is supposed to be applied on hands and body. Apply soap or cleanser on hands or body. Then, gently move the brush on your skin in circles.

SMALL BRUSH



Face brush is made up of very soft bristles, specially designed for tender skin. To achieve the best results of exfoliation, you should first clean your face. Apply desired cleanser to your brush and gently move the brush across your face in circular motion. Please touch the brush on your face in a flat manner instead of placing it in an angle.

PUMICE STONE



Pumice Stone is a light-yet-abrasive stone that aids in keeping skin bright and healthy. It can be used to remove dry and dead skin from the feet with no pain. For effective results, immerse your feet in warm water for a while, then use pumice attachment on your feet by moving it in circles. It can also be used to remove hard cells from elbow or from the foot sole.

FACE SPONGE



Cosmetic soft sponge gently massages your skin to make it radiant and fresh. It works efficiently to remove makeup, dead skin and dirt from your face. It facilitates absorption of your cream or lotions to nourish your skin. It aids in removal of fine lines and black spots. After deeply exfoliating and cleaning your face, dip the sponge in the desired face cream. Gently move the sponge across your face in circular motion.

SETUP GUIDE

Preparing the Unit for Usage

4 pcs AA batteries power the unit. For a longer operating time, it is recommended to use alkaline batteries. Please make sure that your hands are dry while inserting batteries.

Inserting the Battery

1. Use one hand to hold the unit downward and the other hand to hold battery cap with thumb at the "open" marked position.
2. Pull off the cap after pressing for a long time
3. Insert 4 pcs AA batteries.
4. Replace battery cap.
5. Apply different applicators (small brush, sponge, pumice stone and big brush) to the unit for different usage and press the control button (on Main unit) to start the treatment.



Note: Make sure to check the polarity while inserting the batteries.

FAQs

Can I use it for face massage?

Answer: Yes, you can use it for face massage. It is designed for facial massage, scrubbing and cleansing.

How does it work?

Answer: It comes with a battery slot and works with 4 AA batteries.

Is it waterproof?

Answer: No, it is not. But it is water resistant and safe to use in the bathroom.

Can I use it as a body sponge?

Answer: You can also use it on the whole body. For that, just attach the body sponge and rub it in circular motion with soap or body lotion all over the body.

Can I use it as a makeup blender?

Answer: No, it does not work well as a makeup blender. It is a flat cleansing sponge with a rugged surface.

How many times can use this device?

Answer: This is reusable. You can use it every time whenever you clean your face in a day.

Can I use it to remove makeup?

Answer: Yes, you can do that with sponge but ideally, it is meant for use while washing your face.

Does it harm skin?

Answer: No, does not, it even suits well if your skin is sensitive.

BENEFITS OF USING THE FACIAL BRUSH

- It can be used on different body parts to remove the top layer of dead skin.
- Removal of dead skin improves the skin tone.
- It stimulates blood circulation and makes skin firmer and healthier.
- Regular massaging reduces wrinkles.
- It is safer than anti-aging creams as it does not require use of synthetic acids.

USING THE UNIT SAFELY

Maintenance and Replacements

- Clean the unit and all applicators after use, each time.
- Use damp cloth to wipe clean any residues remaining on unit body.
- Rinse applicators with tap water. Max water temperature should be 50°C or 122°F.
- Rinse the brushes thoroughly from inside to outside.
- Let air dry all the parts.

- Store the unit and parts back to the storage box for next usage.
- Replace the brushes every 2 to 3 months for effective cleaning and exfoliation results.
- Replace the brush if the deformed filaments can no longer be restored to their original shape.
- Worn filaments would irritate your skin and tend to spread impurities.
- Do not use any alcohol, acetone, petrol, and abrasive cleaning agents etc, to clean the plastic parts.

WARNINGS

- Do not use it on broken skin or areas of the skin with rash, severe acne or other skin ailments.
- Use a headband to keep your hair away from your face.

Important: If the unit is to be disposed for long time, make sure to take out the batteries.

BOX CONTAINS

This pack contains:

- Massager main unit
- Soft sponge cleaning head
- Big size cleaning brush head
- Pumice stone head
- 4 AA batteries



4 AA batteries



Big Brush



Small Brush



Pumice Stone



Face Sponge

TECHNICAL SPECIFICATIONS



Power Source	4*AA batteries provided in the box
Dimensions	L-8 cm, W -4cm, H -18cm
No. of Applicators	4,(small brush, big brush, sponge, pumice stone surface)

CUSTOMER SUPPORT



CONTACT ADDRESS

USA

NURECA INC.USA

276 5th Avenue, Suite 704-397, New York (NY) - 10001, USA

INDIA

Corporate Office (Mumbai)

Nectar Biopharma Pvt Limited

Quest Offices, The Parinee Crescenzo, C38-39, G Block,

8th floor, B Wing, Bandra Kurla Complex, Bandra(E)

Mumbai, Maharashtra 400051

Call On

USA: +1 212-634-4563

India (Toll free): 18002709565

Website: www.drphys.io

Email: shiatsu@drphys.io

Connect with us on social networks

Facebook: [drphysiousa](#)

Instagram: [drphys.io](#)

Youtube: [NurecaUsa](#)

COPYRIGHT © 2018 DR PHYSIO. ALL RIGHTS RESERVED